




CE approved
ISO 9001 : 2000

CROWN

MERCURY-FREE SPHYGMOMANOMETER

LCD Instead of Mercury
No Pollution
Measurement Accuracy
Suitable for Hospitals, Clinics and Family

CROWN



Two Measuring
1. Auscultatory
2. Oscillography



Measuring Range : 0-300 mmHg
Pulse Range : 30 to 199bpm
Accuracy
— Pressure : ± 3 mmHg
— Pulse : $\pm 5\%$
Battery Type : AA Type
Adapter (Optional) : 4.5v DC/300mA.
Dimensions : 360x96x66 mm
Weight : 700g



CROWN

MERCURY-FREE
SPHYGMOMANOMETER

CE approved
ISO 9001 : 2000



CROWN MERCURY-FREE SPHYGMOMANOMETER

The mercury-free sphygmomanometer opens up a completely new era of professional blood pressure measurement. Recently, awareness of the issue of environmental conservation has been increasing and the use of mercury has become a world-wide concern. Mercury is a toxic metal that can poison whole communities and remain in the local environment for many years, creating serious health risks. Therefore,

the banning of mercury is essential in order to maintain a clean and safe environment.

The CROWN Mercury free B.P. Apparatus was designed to be used as an alternative to the mercury column, while at the same time preserving the local environment and creating an easy-to-use product for hospital personnel.

Reading :

The result will be shown on the screen in black column shape between the SYS and DIA, automatically, so there will be no use of any stethoscope for measurements.

Display : Liquid crystal digital display with SYS, DIA & PULSE display at the same time on the screen. • Power : 2AA Battery • Dimension : 360mm x 96mm x 66mm



MEDI WAVES INC.

H.O. : B-68, G.T. Karnal Road, Industrial Area, Delhi-110033 (India)

Works : Plot No. 92, Sector-57, Phase IV, HSIIDC Kundli,
Sonapat, Haryana-131001 (India)

Helpline No. : 011-42384227, 42384365

E-mail : sales@medi-waves.com Website : www.medi-waves.com

ISO 9001:2015 CERTIFIED
WHO-GMP CERTIFIED
ISO 14001:2004 CERTIFIED
OHSAS 18001:2007
ISO 13485 : 2016 CERTIFIED
CE CERTIFIED