



Model No. : P-010

Oxygen Hood

A device consisting of a rigid/semi-rigid transparent plastic shell that forms an enclosure over an infant's whole body, or the head only, in order to provide an enriched environment of oxygen (O₂) to increase the patient's O₂ uptake. It is connected to an O₂ source and may be used concurrently with increased humidification and temperature control. It is designed to be used for patients adverse to oxygen delivery devices such as a nasal cannula or face mask. This device may include the tubing, a diffuser (to disperse the flow of incoming O₂), O₂ concentration This is a reusable device.



SMALL



MEDIUM



LARGE

FEATURES

- Single Piece Round shape
- Made up of autoclavable Imported Polycarbonate material.
- Trauma free Silicon neck adjustable flap.
- Bi-lateral Oxygen Nozzle, Prevents direct flow of cold Oxygen on Patients head.
- Sizes Available : Small, Medium & Large

Technical Specifications



MEDI WAVES INC.

ISO 9001:2015 CERTIFIED

WHO-GMP CERTIFIED | ISO 14001:2004 CERTIFIED | OHSAS 18001:2007

ISO 13485 : 2016 CERTIFIED | CE CERTIFIED

H.O. : B-68, G.T. Karnal Road, Industrial Area, Delhi-110033 (India)

Works : Plot No. 92, Sector-57, Phase IV, HSIIDC Kundli,
Sonapat, Haryana-131001 (India)

Helpline No. : 011-42384227

E-mail : sales@medi-waves.com

Website : www.medi-waves.com

ISO 9001:2015 CERTIFIED

WHO-GMP CERTIFIED | ISO 14001:2004 CERTIFIED | OHSAS 18001:2007

ISO 13485 : 2016 CERTIFIED | CE CERTIFIED